



SPRINGFIELD CRIB



Our updated Springfield Crib offers the functionality of our Standard Crib, with additional protection for an active child.

This model allows for the sides and head end to drop for access to the patient by multiple care givers, along with an add on top with vinyl safety panels, providing extra safety for the active patient.

The height of the updated crib rails are slightly taller, while the vinyl panels are shorter, keeping the overall crib clearance the same.


FEATURES INCLUDED

- Manual Head Elevation to 35 Degrees
- Manual Knee Elevation to 15 Degrees
- 13 cm Pressure Reduction Mattress
- 2 Drop Sides & Drop Head End Rail
- Add-On Crib Top - (Shown Bottom Right)
Note: Vinyl Ceiling Panel Optional
- IV Pole
- Shelf with Oxygen Tank Cut Outs
- 13 cm Heavy Duty Casters
(2 Steer, 2 Brake)
- Teething Bars
- Wall Bumper System
- Confetti Epoxy Finish
- Steel Construction & Factory Assembly


**Recommended for Patients
up to 68 kg Maximum**

OPTIONS & ACCESSORIES

- Vinyl Ceiling Panel
- Additional IV Poles - Up to 3
- Portland or Dialysis Access
- Trendelenburg Function
- Stainless Steel
- Oxygen Tank Holder
- Fitted Sheets
- Additional Epoxy Colors Available



**Ceiling
Panel -
NOT
Included**



**Add-On
Crib Top -
Included**

PC1935-T

Crib Size No Bumpers: 76x137 cm
Crib Size With Bumpers: 89x150 cm
Mattress Size: 71x132 cm
Overall Crib Height: 190 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Height of Add-On Crib Top: 43 cm
Crib Weight: 107 kg

PC1943-T

Crib Size No Bumpers: 76x152 cm
Crib Size With Bumpers: 89x165 cm
Mattress Size: 71x147 cm
Overall Crib Height: 190 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Height of Add-On Crib Top: 43 cm
Crib Weight: 116 kg

PC1955-T

Crib Size No Bumpers: 91x152 cm
Crib Size With Bumpers: 104x165 cm
Mattress Size: 86x147 cm
Overall Crib Height: 190 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Height of Add-On Crib Top: 43 cm
Crib Weight: 124 kg

PC1975-T

Crib Size No Bumpers: 91x183 cm
Crib Size With Bumpers: 104x196 cm
Mattress Size: 86x178 cm
Overall Crib Height: 190 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Height of Add-On Crib Top: 43 cm
Crib Weight: 136 kg

AVAILABLE COLORS



ISO 13485:2016
Registered by
EAGLE Registrations Inc.

