



STANDARD CRIB



Our Standard Crib provides long-lasting, high-quality construction and all the basic functionality needed to be cost effective and versatile for general use with newborns through young children in the hospital.

Crib tops can be added onto the Standard Crib after initial purchase, providing additional flexibility if the needs change within the hospital.

FEATURES INCLUDED

- Manual Head Elevation to 35 Degrees
- Manual Knee Elevation to 15 Degrees
- 13 cm Pressure Reduction Mattress
- 2 Drop Sides
- IV Pole
- Shelf with Oxygen Tank Cut Outs
- 13 cm Heavy Duty Casters (2 Steer, 2 Brake)
- Teething Bars
- Wall Bumper System
- Confetti Epoxy Finish
- Steel Construction
- Factory Assembly

**Recommended for Patients
up to 68 kg Maximum**

OPTIONS & ACCESSORIES

- Add-On Crib Top - Shown Below
- Vinyl Ceiling Panel - For Add-On Crib Top
- Additional IV Poles- Up to 3
- Portland & Dialysis Access
- Trendelenburg Feature
- Stainless Steel
- Oxygen Tank Holder
- Fitted Sheets
- Additional Epoxy Colors Available



PC1913

Crib Size No Bumpers: 76x112 cm
Crib Size With Bumpers: 89x124 cm
Mattress Size: 71x107 cm
Overall Crib Height: 147 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Crib Weight: 77 kg

PC1931

Crib Size No Bumpers: 76x137 cm
Crib Size With Bumpers: 89x150 cm
Mattress Size: 71x132 cm
Overall Crib Height: 147 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Crib Weight: 89 kg

PC1942

Crib Size No Bumpers: 76x152 cm
Crib Size With Bumpers: 89x165 cm
Mattress Size: 71x147 cm
Overall Crib Height: 147 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Crib Weight: 96 kg

PC1952

Crib Size No Bumpers: 91x152 cm
Crib Size With Bumpers: 104x165 cm
Mattress Size: 86x147 cm
Overall Crib Height: 147 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Crib Weight: 112 kg

PC1973

Crib Size No Bumpers: 91x183 cm
Crib Size With Bumpers: 104x196 cm
Mattress Size: 86x178 cm
Overall Crib Height: 147 cm
Floor to Top of Mattress
Platform: 79 cm
Height of Side Rails: 66 cm
Crib Weight: 129 kg

AVAILABLE COLORS



ISO 13485:2016
Registered by
EAGLE Registrations Inc.

